

MESI mTABLET Thermometer

Non-Contact and Contact Thermometer in One

Perform contactless and diverse contact temperature measurements without manual data entry.

VIEW, SAVE OR DISCARD RESULTS ON MEDICAL-GRADE TABLET

- > Fast and easy-to-use **contactless forehead option**
- > **Contact thermometer** for specialised measurements and/or confirmation
- > **Advanced analysis options** with patient history

WIRELESS & PORTABLE



DISTANCE ALIGNMENT TOOL



ARMPIT, RECTAL OR ORAL MEASUREMENTS



SAVE DIRECTLY INTO EHR



SHARE FOR A SECOND OPINION



Why MESI mTABLET Thermometer?

Non-contact forehead thermometer

- Measured temperature from 34.0 °C to 42.0 °C.
- Automatic distance alignment tool.
- Colour result indicator for immediate feedback.
- Minimum cross-contamination risk.
- Smoother and less invasive process with less distress for the patient.
- No disposable probe covers needed – enhanced sustainability and lower total cost of ownership.



Contact thermometer

- Measured temperature from 34.0 °C to 42.0 °C.
- Armpit, oral or rectal measurements.
- Disposable contact probe shields for reduced cross-contamination risk.



Application extensions

The **MESI mTABLET Thermometer** is not only an advanced temperature-measuring device. You can also use it within the MESI mTABLET Vitals (with heart rate, blood pressure and SpO2) for complete assessment of vital signs. By adding smart apps like Protocol, you can also make temperature a part of custom protocols at your practice.

TEMP

Temperature

Vitals

Vital signs

Protocol

Health assessment protocol

Photo

Camera

Worklist

Patient worklist



MESI, Ltd.
Leskoškova cesta 11a
1000 Ljubljana,
Slovenia, EU

E: info@mesimedical.com
T: +386 1 620 34 87

www.mesimedical.com

f MESIdoo

🐦 MESImedical

in MESI



European production
and development



Worldwide presence

FDA

Food and Drug
Administration
cleared

MDR

EU Medical Device
Regulations
compliant



ISO 9001 Q-1664
ISO 13485 M-049

ISO 9001 and
ISO 13485 certified



MDSAP compliant